

Modulator

EMERGENCY
RESPONSE BACKPACK



Modulator Backpack
has won the
**MPA Health &
Safety Award** for
Safety Operations
On-Site



The Ultimate **Emergency Response Backpack**

Unlock a new era of safety and preparedness with the Modulator Emergency Response Backpack meticulously crafted for professionals in the rail, wind farm, mining, heavy machinery, and oil & gas industries. Engineered to perfection, it goes beyond expectations, ensuring you're ready to respond with confidence in any emergency.



Thermal Insulation



Foul Weather Proof



Working Load Limit



Lightweight Design



Quick Release Buckle



Flexible storage for either a **ViVest** or **HeartSine** defibrillator

- ✔ ViVest or HeartSine Defibrillator
- ✔ Burns Module
- ✔ Eye Module
- ✔ Trauma Module
- ✔ Dressings and Bandage Module
- ✔ Cuts and Grazes Module



Fabric complies with **RIS-3279-TOM** High-Visibility Standard



Includes a **Personal First Aid Kit** in the quick access pocket



“What first aid provision is inside the backpack?”



With a ViVest / HeartSine AED*

Can save someone from a sudden cardiac arrest incident. *Optionally available either with a X-Series or Samaritan defibrillator.



CONTENTS	QTY
AEROBURN STERILE FACE GEL DRESSING FACE MASK	1
AEROBURN STERILE GEL DRESSING 20x20CM	1
AEROPLAST EMERGENCY FOIL BLANKET SILVER 140x210CM	1
AEROGLOVE NITRILE GLOVES POWDER FREE MEDIUM – PAIR	1
AEROGLOVE NITRILE GLOVES POWDER FREE EXTRA LARGE – PAIR	1



BURNS MODULE



TRAUMA MODULE

CONTENTS	QTY
AEROWOUND TRAUMA DRESSING MEDIUM 10x18CM PAD	1
AEROBAND CC COTTON CALICO TRI-BANDAGE 90x90x127CM	1
AEROPLAST EMERGENCY FOIL BLANKET SILVER 140x210CM	1
AEROGLOVE NITRILE GLOVES POWDER FREE LARGE – PAIR	1
RAPIDSTOP TOURNIQUET – ORANGE STRAP	1
AEROWOUND TRAUMA DRESSING LARGE 15x18CM WOUND PAD	1

Includes a RapidStop® Tourniquet

Relies on key technologies to achieve rapid, intuitive and easy one-hand self-application.



First Aid Provision continues on following page



CUTS AND GRAZES MODULE

CONTENTS	QTY
AEROPLAST WASHPROOF ASSORTED PLASTERS – ENVELOPE x 10	1
AEROWOUND TRAUMA DRESSING MEDIUM 10x18CM WOUND PAD	1
AEROWOUND TRAUMA DRESSING LARGE 15x18CM WOUND PAD	1
AEROGLOVE NITRILE GLOVES POWDER FREE MEDIUM – PAIR	1
AEROGLOVE NITRILE GLOVES POWDER FREE LARGE – PAIR	1
AEROPORE CLEAR ADHESIVE TAPE 2.5CMx5M	1
AEROHAZARD CLINICAL WASTE BAG 381x711x990MM – PACK x 25	1
AEROWIPE SALINE CLEANSING WIPES – EACH	10

CONTENTS	QTY
AEROWASH EYEWASH 250ML – EACH	1
AEROBAND CC COTTON CALICO TRI-BANDAGE 90x90x127CM	1
AEROGLOVE NITRILE GLOVES POWDER FREE MEDIUM – PAIR	1
AEROGLOVE NITRILE GLOVES POWDER FREE EXTRA LARGE – PAIR	1
AEROWOUND HSE EYE PAD – EACH	3
AEROWASH 20ML EYEWASH – EACH	2



EYE MODULE



DRESSINGS AND BANDAGE MODULE

CONTENTS	QTY
AEROWOUND TRAUMA DRESSING LARGE 15x18CM WOUND PAD	1
AEROBAND CC COTTON CALICO TRI-BANDAGE 90x90x127CM	1
AEROGLOVE NITRILE GLOVES POWDER FREE MEDIUM – PAIR	1
AEROGLOVE NITRILE GLOVES POWDER FREE LARGE – PAIR	1
AEROWOUND TRAUMA DRESSING MEDIUM 10x18CM WOUND PAD	1
AEROWOUND HSE MEDIUM DRESSING 12x12CM	1
AEROWOUND HSE LARGE DRESSING 18x18CM	1

CONTENTS	QTY
AERO FIRST AID AT WORK GUIDANCE LEAFLET	1
AEROBAND NON-WOVEN TRI-BANDAGE 90x90x127CM	2
A-CARE WASHPROOF ASSORTED PLASTERS – ENVELOPE x 10	1
AEROPIN SAFETY PINS 3 SIZES NICKLE PLATED – PACK x 6	1
AEROWIPE SALINE CLEANSING WIPES – EACH	2
AEROWOUND HSE LARGE DRESSING 18x18CM	1
AEROGLOVE VINYL GLOVES POWDER FREE LARGE – PAIR	1



PERSONAL FIRST AID KIT



INCLUDING A USER-FRIENDLY DEFIBRILLATOR

ViVest PowerBeat AEDs


Semi-automatic defibrillators (X3 or X1) designed for anyone to operate, offering quick and safe emergency response through visual guidance, voice prompts, and pre-connected smart pads.



With ViVest X3 AED feat. 5" LCD Display

- ✔ Animated step-by-step visual instructions
- ✔ One-button activation plus CPR metronome
- ✔ 5 year battery life (up to 200 shocks)
- ✔ Bluetooth and 4G remote connectivity
- ✔ Automated ECG recording (up to 8 hours)
- ✔ Daily automated self-testing
- ✔ Single Paediatric mode button
- ✔ IP55 dust and water resistance
- ✔ Supports up to 3 selectable languages

Aero Code: AK-MOD-RS-X3



With ViVest X1 AED feat. Pictograms

- ✔ A budget-friendly model offering all the features of the X3, but with a colour pictogram visual panel instead of an LCD display.

Aero Code: AK-MOD-RS-X1



HeartSine Samaritan AEDs

A range of semi-automatic (500P or 350P) or fully automatic (360P) defibrillators designed to deliver reliable emergency care with user-friendly guidance and a compact design.



With HeartSine 500P feat. CPR Advisor

- ✔ CPR Advisor: real-time CPR feedback
- ✔ Semi-Automatic Shock Delivery
- ✔ Easy to Follow Visual and Verbal Prompts
- ✔ IP56 Rated Ingress Protection
- ✔ Advanced Technology
- ✔ Pad-Pak: Battery & Electrodes combined
- ✔ Portable and Lightweight at just 1.1 kg
- ✔ Software updates via USB
- ✔ Free replacement Pad-Pak initiative

Aero Code: **AK-MOD-RS-500**



With HeartSine 350P Semi-Automatic AED

- ✔ A budget-friendly semi-automatic AED with all the features of the 500P except for CPR Advisor. **This defibrillator will prompt the user to administer an electrical shock if required.**

Aero Code: **AK-MOD-RS-350**



With HeartSine 360P Fully Automatic AED

- ✔ A budget-friendly fully automatic AED with all the features of the 500P except for CPR Advisor. **This defibrillator will automatically deliver an electrical shock if needed.**

Aero Code: **AK-MOD-RS-360**