





## FABRIC PLASTERS

✓ MDD Compliant



Designed to treat a variety of injuries, **A-Care Economy Latex-Free Fabric Dressings** provide an effective solution for all types of wounds. Manufactured from light, elasticated fabric, each latex-free dressing is strong and secure and can withstand regular handwashing to protect the injury. The dressings stretch, bend, flex, and conform to body contours for maximum adhesion, especially over awkward areas such as fingertips and knuckles. This prevents lifting and edge curling and keeps the wound clean and comfortable. These durable fabric plasters are both hypo-allergenic and latex free and can be used by a diverse range of patients including those who suffer from skin complaints or allergies.

### AVAILABLE IN:



ACF100 – Strip (Box 100)



ACF101 – Large Patch (Box 50)



ACF102 – Fingertip (Box 50)



ACF105 – Small Patch (Box 100)



ACF106 – Assorted (Box 100)

