

AERO WIPE™

CLEANSING SWABS

ALCOHOL

- ✓ MDD Compliant
- ✓ UKCA Marked & CE Marked

AW7100_UK



AeroWipe™ Skin prep wipes contain 70% Isopropyl Alcohol which is an effective antibacterial formula for cleansing the skin before injection. Alcohol acts as a disinfectant to clean and effectively remove bacteria, perfect for preparing and cleansing skin for procedures such as blood collection or inoculation of a vaccine. AeroWipe™ are UKCA and CE-marked, manufactured in accordance with ISO 13485 quality standards. Each wipe is 3cm x 3cm in size, and individually wrapped in a sealed sachet, ensuring the maintenance of sterility before use.

AVAILABLE IN:



Alcohol



3cm x 3cm (Box 100)

Alcohol



3cm x 3cm