



AERO TACTICAL™

TRAUMA MEDIC BERGEN



Modular &
Customisable



Optimised Visibility
& Accessibility



Built for
Endurance



Waterproof



Fits larger
CO₂ cylinder

DESIGNED WITH FORCE EXPERTISE, FOR OPERATIONAL EXCELLENCE

Elevate your state of readiness with the AeroTactical™ Trauma Medic Bergen – a meticulously crafted solution, engineered through extensive trials with various police forces across the UK. The result of this intensive product development phase is a customisable and exceptionally robust bergen, capable of accommodating modular and configurable internals. Engineered for swift frontline deployment, the AeroTactical™ Trauma Medic Bergen empowers confident response in every critical moment.

 **Factor of Safety 7:1**
conforms to EN1492-1

MEDIC BERGEN: KEY BENEFITS

- ⚙️ Designed by forces for forces.
- ⚙️ Modular and configurable inserts.
- ⚙️ Front pouch fits any defibrillator model →
- ⚙️ Hard-wearing (reinforced base and sides).
- ⚙️ Freestanding and upright on bottom or side.
- ⚙️ Zips can be security sealed.
- ⚙️ Chest harness increases carrying comfort.



 **Working Load Limit: 30 kg**

- ❖ Low light visibility with reflective strips →
- ❖ Waterproof, providing ingress protection.
- ❖ External surfaces easily wipe clean.
- ❖ Able to fit larger CO² cylinder.
- ❖ High visibility fleck on zipper pull cords.
- ❖ Document sleeve in the back panel.
- ❖ Available with 3 different coloured straps.



❖ Reinforced rubber base

❖ VELCRO® attachment straps allows quick release of modular inserts



❖ Reassuringly waterproof



EMS-PAN-BAG-BLACK

AeroTactical™ Trauma Medic Bergen – BLACK

AVAILABLE WITH DIFFERENT COLOURED STRAPS



EMS-PAN-BAG-GREEN

AeroTactical™ Trauma Medic Bergen – GREEN



EMS-PAN-BAG-RED

AeroTactical™ Trauma Medic Bergen – RED

MEDIC BERGEN INSERTS: KEY BENEFITS

- ✓ Ability to insert ID tags onto pouches.
- ✓ Pouches have clear display panel to quickly identify contents.
- ✓ Quick release zips with pull cords.
- ✓ Pouch and insert surfaces easily wipe clean.
- ✓ High visibility lining to quickly identify any missing pouches, etc.
- ✓ VELCRO® backing allows quick release of pouch / insert.
- ✓ Waist buckle available for added comfort and security.
- ✓ First aid kit contents can be custom packed by Aero.



EMS-PAN-IP-L-BLUE

AeroTactical™ Enclosed Zippered
Inner Pouch – Large BLUE



EMS-PAN-IP-L-RED

AeroTactical™ Enclosed Zippered
Inner Pouch – Large RED



EMS-PAN-IP-O2

AeroTactical™ O2/N2O3
Cylinder Sleeve




EMS-PAN-RBI

AeroTactical™ Removeable
Bungee Insert



EMS-PAN-IP-L-BLACK

AeroTactical™ Enclosed Zippered
Inner Pouch – Large BLACK



EMS-PAN-ELI

AeroTactical™ Removeable
Elastic Loop Insert



EMS-PAN-IP-M

AeroTactical™ Enclosed Zippered
Inner Pouch – MEDIUM



EMS-PAN-IP-S


AeroTactical™ Enclosed Zippered
Inner Pouch – SMALL

MEDIC BERGEN RANGE: TECHNICAL SPECIFICATIONS



AeroTactical™ Trauma Medic Bergen

Height: **59 cm**
 Width: **40 cm**
 Depth: **28 cm**
 Internal Volume: **c. 46 L**
 Working Load Limit: **30 kg**
 Factor of Safety **7:1**
 Conforms to **EN1492-1**



Outer Main Material:
840D Nylon with TPU coating
 Interior Lining Material:
600D Polyester with PU coating
 Padding (different size foams are used):
EPE foam 3mm/ 5mm/8mm; EVA foam 3mm/4mm; PE board 1mm/2mm




Webbing Material: **Polyester**
 Zipper: **SBS double-runner, quick release zips with pull cords**

AeroTactical™ Enclosed Zippered Inner Pouch – Large



Height: **38 cm**
 Width: **18.5 cm**
 Depth: **8 cm**
 Material:
500D PVC Tarpaulin

AeroTactical™ Removeable Elastic Loop Insert




H: **33 cm**
 W: **14 cm**
 D: **0.5 cm**
 Material:
500D PVC Tarpaulin, PE board

AeroTactical™ O2/N2O3 Cylinder Sleeve



H: **38 cm**
 W/D: **12.5 cm**
 Material:
PVC outer, 210D with PU coated lining, EPE foam

AeroTactical™ Removeable Bungee Insert



H: **33 cm**
 W: **14 cm**
 D: **1.5 cm**
 Material:
500D PVC Tarpaulin, PE board

AeroTactical™ Enclosed Zippered Inner Pouch – Medium



H: **18.5 cm**
 W: **18.5 cm**
 D: **8 cm**
 Material:
500D PVC Tarpaulin

AeroTactical™ Enclosed Zippered Inner Pouch – Small



H: **12 cm**
 W: **18.5 cm**
 D: **8 cm**
 Material:
500D PVC Tarpaulin

Contact Aero Healthcare today

For further information on our first aid and trauma provision, defibrillators and how you could benefit, simply get in touch with a member of our team who will be happy to assist you.



HEALTHCARE | WE MAKE IT BETTER

CALL **01403 790 704**

EMAIL **sales@aerohealthcare.co.uk**

VISIT **aerohealthcare.com/uk**

© Aero Healthcare