




BEURER LIFE PAD® RESUSCITATION AID

- ✓ MDR Compliant
- ✓ UKCA Marked & CE Marked



Order code: **LIFEPAD**

A cardiac arrest can affect anyone, at any time. The award-winning **Beurer LifePad®** supports laypersons step-by-step, ensuring safe and effective CPR (cardiopulmonary resuscitation). It's designed for ease of use with no prior training required. Red markings show proper chest placement, while a raised surface guides hand positioning. Made of flexible material to fit body contours, the LifePad® is suitable for users aged 12 and older. Its acoustic and visual feedback system provides guidance on compression depth and frequency, with a rhythm signal set at 100 beats per minute – compress with each beep and release after. Colour-coded LEDs indicate correct compression pressure, illuminating fully when sufficient pressure is applied. This combination of visual and acoustic feedback offers clear guidance and reassurance, even in high-stress emergency situations, helping anyone perform life-saving CPR measures effectively. Can also be used in combination with a defibrillator (audible signal tone can be switched off for this purpose).



AERO
Healthcare