



## CR Plus Charge-Pak Charging Stick With Electrodes

- ✓ MDD Compliant
- ✓ UKCA Marked & CE Marked



The Physio-Control Lifepak CR Plus Replacement Kit is designed to ensure the continued reliability and functionality of your defibrillator. With a stand-by life of 2 years and a replace-by date aligned with this period, the kit includes essential components for seamless replacement and maintenance. It features a CHARGE-PAK™ Battery Charging Unit to keep your defibrillator's battery fully charged, a set of Quickpak Electrodes (available with either 1 or 2 pairs), and a discharger for the safe disposal of used CHARGE-PAK batteries. The Quickpak Electrodes, specifically engineered for use with the Physio-Control Lifepak CR Plus Fully Auto and Semi-Auto Defibrillators, have a 2-year lifespan to ensure optimal performance when needed. This comprehensive kit helps you maintain readiness and compliance, providing peace of mind with the assurance that all components are current and functional.

### AVAILABLE IN:



11403-000002



One Electrode

11403-000001



Two Electrodes

For this Defibrillator:



Lifepak CR Plus