

COHESIVE BANDAGE

- ✓ MDD Compliant
- ✓ UKCA Marked & CE Marked



The AeroBan™ Cohesive Bandage affords extra comfort as it doesn't pull or attach itself to body hair. This bandage offers more support for wearers as it holds itself in place without the need for retention clips, safety pins, or microporous tape. AeroBan™ Cohesive Bandages are premium quality and comfortable to use. The bandage can be easily manipulated to the size required and can adapt to body contours. Its self-adhesive design allows it to stick to itself without sticking to the skin. The bandage is designed for supporting joints without discomfort. Additionally, the AeroBan™ Cohesive Bandage is also latex free, which means it can be used by individuals that suffer with certain skin complaints, such as allergies. It is classed as non-sterile.

AVAILABLE IN:



5cm x 4m



7.5cm x 4m