

ELASTIC TUBULAR BANDAGE

- ✓ MDD Compliant
- ✓ UKCA Marked & CE Marked



AeroForm™ Tubular Bandages are elasticised support bandages that are ideal for flexible compression and support. Bandages can be cut to desired length or doubled over for extra compression offering firm but flexible support on limbs. AeroForm™ Tubular Bandages are held in place without tying or tapes, resulting in a comfortable long wearing bandage. They provide even support across damaged joints and muscles due to its elastic self supporting properties resulting in a comfortable support bandage. A non-sterile, tubular, elastic fabric (e.g., polyamide, Lycra®) sleeve intended to be externally worn on a limb to provide support or local pressure to a part of the body, especially a joint, for various preventative/therapeutic applications (e.g., support soft tissue injuries) while enabling movement. Sometimes referred to as a tubigrip, it may be designed to accommodate one or more anatomies (e.g., ankles, knees, wrists). This is a single-patient, reusable device.

AVAILABLE IN:



Small Limbs (6.5cm x 1m)



Adult Limbs (7cm x 1m)



Large Adult Limbs (7.5cm x 1m)



Adult Legs (8.5cm x 1m)



Large Adult Legs (10cm x 1m)

