

Product Code: HS-WC-014

# ORANGE HARD CASE FOR HEARTSINE DEFIBRILLATORS




## WEATHERPROOF, LIFESAVING TECHNOLOGY



Ensuring your defibrillator is protected against extreme weather, whilst remote working

- ✓ Tested to IP67
- ✓ Tested to MIL-STD-810F for Vibration and Impact Resistance
- ✓ Crushproof
- ✓ Dustproof
- ✓ Waterproof
- ✓ Lightweight
- ✓ Bespoke Thermal Protection Form



## EXTERNAL DIMENSIONS

Interior (mm): 411  
Width (mm): 322  
Depth (mm): 168



# TECHNICAL SPECIFICATIONS

All hard cases containers meet the applicable requirements of MIL-STD-810F , FED-STD- 101C and ATA-300 (Category 1) as follows:

## BUOYANCY

Floats in salt water with 13.61 kg load

## APPROVALS

IP67 – Defstan 81-41/Stanag 4280 MIL C-4150-J

## PACKAGING

Single Unit 420 x 358 x 182mm, 3.48 kg (with foam),  
3.06 kg (without foam) 3.48 (with dividers)

Master pack of 4 units: 752 x 430 x 374mm, 13.42 kg (with foam),  
10.28 kg (without foam), 14.92 kg (with dividers)

Maximum number on a pallet = 60 units

## ACCESSORIES

3 piece foam set (1 pick n pluck section)

Panel Frame

Divider Set

TSA Lock

Nameplate

## OTHER INFORMATION

Single handle with non-slip overmoulded grip

Easy-open latches


Automatic pressure equalisation valve

Rubber 'O-ring' seal

Solid wall design

Scratch resistant

## ADVANCED PROPRIETARY TECHNOLOGY



The foam we used is a closed cell, cross-linked polyethylene foam manufactured using Zotefoams unique production process with a operating temperature reaching +100 cercius to a miniumum working temperatures of -70.



## GUARANTEE

You break it, we replace it. Forever.

