



IPAD SP1 Defibrillator Battery

- ✓ MDD Compliant
- ✓ UKCA Marked & CE Marked



The battery pack is a crucial component of the defibrillator machine, as an adequate power supply is essential for providing life-saving treatment during a cardiac arrest. The battery has a 5-year lifespan and, when in standby mode at operable temperatures, performs daily, weekly, and monthly self-tests to ensure the device is always rescue ready. Designed exclusively for use with the I-PAD SP1 and SP2 defibrillator, the battery provides a minimum of 200 shocks (150 joules) and features a simple click-and-replace procedure for easy installation. After being clipped to the back of the defibrillator and set to standby mode, the battery initiates regular self-checks of the subsystems, hardware, software, and electrode pads to ensure all components are functioning correctly. This self-check process allows the defibrillator to remain unattended for long periods while maintaining its readiness for emergency use.

CUA1007B



For this Defibrillator:



IPAD SP1
& SP2