



## CR2, Fully Automatic, USB, Handle

- ✓ MDD Compliant
- ✓ UKCA Marked & CE Marked



The Lifepak® CR2 is a fully automatic defibrillator designed for quick, life-saving response during sudden cardiac arrest. It replaces the Lifepak® CR Plus, offering faster heart rhythm analysis and automatic shock delivery when needed. The device includes clear visual instructions, CPR coaching with a metronome for timing, and a convenient handle for easy transport. With a 4-year battery and electrode pad life, it requires minimal maintenance. Additionally, the CR2 allows patient and event data transfer via USB for easy post-event analysis and record-keeping.

PHE1044



- ✓ Easily Transportable
- ✓ Fast Heart Rhythm analysis
- ✓ Simple battery replacement process