



CR2, Fully Automatic, Wi-Fi, Handle

- ✓ MDD Compliant
- ✓ UKCA Marked & CE Marked



The Lifepak® CR2 is a fully automatic defibrillator that replaces the Lifepak® CR Plus, designed to support first responders during sudden cardiac arrest. It quickly analyzes heart rhythms and delivers shocks automatically, prompting bystanders to stand clear for safe use. Known for its rapid response, the CR2 features user-friendly pictorial instructions and CPR coaching with a metronome for proper compression timing, all housed in a transport-friendly design. Equipped with WiFi capabilities, it wirelessly transfers patient and event data for analysis and performs self-monitoring checks to notify managers of maintenance needs. With a 4-year battery and pad life, the Lifepak® CR2 ensures reliability and readiness for emergencies.

PHE1039



- ✓ Easily Transportable
- ✓ Fast Heart Rhythm analysis
- ✓ Simple battery replacement process