



CR2, Semi-Automatic, USB, Handle

- ✓ MDD Compliant
- ✓ UKCA Marked & CE Marked



The Lifepak® CR2 is a semi-automatic defibrillator designed to replace the CR Plus unit, providing first responders with fast and reliable support during sudden cardiac arrest. It quickly analyzes heart rhythms, advising responders when to deliver a shock via a button press. The device includes user-friendly pictorial and visual instructions, CPR coaching, and a metronome to maintain compression timing. A convenient handle ensures easy transport, while its 4-year battery and pad life minimize maintenance. Additionally, event data can be transferred via USB for post-incident analysis, enhancing both usability and performance.

PHE1043



- ✓ Easily Transportable
- ✓ Fast Heart Rhythm analysis
- ✓ Simple battery replacement process