

CARDIACT™ LOCKED POLYCARBONATE EXTERNAL CABINET

- ✓ UKCA Marked
- ✓ CE Marked



Shield your defibrillator from the elements using the CardiAct™ Locked Polycarbonate External Cabinet for added protection.

- IP66 rating, sealed against high pressure water hoses and dust
- Suitable for coastal locations
- Fully corrosion resistant in UV stabilised high impact polycarbonate, the same material used in Police riot shields. The UV resistance protects against the harmful effects of the sun. It is also rated to UL94V0 which means they will not support a fire and are self-extinguishing, do not give off toxic fumes.
- Marine grade stainless steel keypad lock, salt spray tested for 10,000 hours
- Internal LED motion sensing light which can be removed and used as a torch during deployment
- Thermostatically controlled heater ensures the defibrillator is kept in optimum conditions
- Internal hook accommodates any defibrillator without a carry case
- Large viewing window allows easy checking on the status of the AED defibrillator
- Internal hinges that do not allow for external tampering
- Double skin clamshell design provides additional insulation in cold weathers
- Comes complete with comprehensive installation instructions and a fixing kit
- Cabinet dimensions 465mm x 435mm x 240mm
- 10-year warranty on the shell of the cabinet; 5-year warranty on the lock against any mechanical malfunction.