

NASOPHARENGEAL AIRWAY

- ✓ MDD Compliant
- ✓ UKCA & CE Marked



The **Nasopharyngeal Airway** is a thin, clear, flexible tube that is inserted into a patient's nostril. The purpose of the NPA is to bypass upper airway obstruction at the level of the nose, nasopharynx or base of the tongue. It also prevents the tongue falling backward on the pharyngeal wall to prevent obstruction. Stays-in-situ without the need for ribbons or safety pins. The Airway Nasopharyngeal has a unique rectangular flange that sits comfortably for the patient and ensures the device cannot be inhaled.

- One-piece design, no parts to assemble
- Designed with a unique and innovative rectangular flange
- Soft material for ease of insertion
- Made of clear material to enable blockages to be identified

AVAILABLE IN:



Size 6



Size 7



Size 8