

NON-WOVEN TRIANGULAR BANDAGE

- ✓ MDD Compliant
- ✓ UKCA Marked & CE Marked



Essential for any first aid kit, our AeroBand™ Non-Woven Triangular Bandage offers support to first aiders. This non-sterile, non-woven triangular bandage can be used as an immobilisation tool, whilst also ensuring support for arms when used as a sling. The function of this triangular bandage can also be used for shoulders and ribs. This bandage is durable and can offer significant support to a casualty whilst expert medical care is on its way. It can also be used as a large dressing to keep sterile pads in place or to fully support injured or wounded limbs. Designed for single use, the AeroBand® Non-Woven Triangular Bandage must be disposed of after it has been used on a casualty, even if it has only been used as a sling. Suitable for immobilisation and limb support.

AVAILABLE IN:



90cm x 90cm x 127cm



CC Cotton Calico Tri-Bandage
90cm x 90cm x 127cm