



AERO PLAST™

PLASTIC

AeroPlast™ Plastic offer an economical, washproof range of conforming bandages perfect for protecting wounds from dirt and germs that can cause infection.



PROTECTIVE WRAPPER

A cold-sealed paper wrapper ensures a sterile bandage, with distinct opening tabs for easy fumble-free opening.

Individually wrapped for guaranteed quality and tamper protection.

FEATURES


- Washproof
- Conforming
- Sterile
- Latex Free
- Hypoallergenic
- SupaBond™ Adhesive Technology




PLASTIC PLASTERS

PRODUCT INDICATIONS

The AeroPlast™ Plastic range is available with a selection of shapes and sizes to ensure the perfect fit in different applications.



**THE ECONOMICAL, WASHPROOF, SHEER AESTHETIC DRESSING
RECOMMENDED FOR ALL GENERAL APPLICATIONS.**



AEROPLAST™ PLASTIC LAYER

A unique plastic layer that bends, stretches and conforms to different shapes and applications. A washproof layer that protects the wound from dirt and germs, while being air permeable allowing it to breathe.

SUPABOND® ADHESIVE

Exclusively engineered for unrivaled latex free and hypoallergenic adhesion.

WOUND PAD

Super absorbent non-woven wound pad that helps reduce skin maceration & aid in faster healing.

NON STICK BARRIER

A non-adherent barrier, provides wound protection and minimises trauma on removal.

PROTECTIVE FILM

PE coated & silicone release papers ensure easy removal without touching the wound pad & helps maintain superior adhesion.