

5th Generation Product Data Sheet

7L

Our Sharpsafe containers are specifically designed to enable the safe disposal of clinical sharps and minimise the risk of associated injuries.

Our 7L is designed for low to medium use areas and ideal for acute and community settings.

Available with different colour lids for easily identifiable waste segregation. Products come with a range of compatible accessories for increased stability, security and safety.

Sharpsafe's new 5th Generation containers launched in 2022 and were the world's first full range of sustainable single use sharps containers. Engineered using low carbon footprint recycled material showing our commitment to championing change.



SUPPORTING
A CIRCULAR
ECONOMY



REDUCING
YOUR CARBON
FOOTPRINT



REDUCING VIRGIN
PLASTIC ENTERING
THE SUPPLY CHAIN



MADE USING LOCALLY
SOURCED LOW
CARBON FOOTPRINT
RECYCLED MATERIALS



Product Certification	
ISO 23907-1:2019	✓
UN3291	✓
BSI Kitemark	✓

Site Certification	
ISO 9001 - Quality Management	✓
ISO 14001 - Environmental Management	✓
ISO 45001 - Occupational Health & Safety	✓

Product Details	
Size:	7L
Usable Capacity:	5.78L
Usable Volume:	85%
Weight:	401g/410g
Fill Line Height:	308mm
Autoclavable:	Yes
Autoclavable Conditions:	134°C
Product Dimensions (mm):	380 x 175 x 175 (HxWxD)
Aperture Dimensions (mm):	110 x 60 (LxW)
GTIN:	5056271506173
Lot/Batch Number:	Not Applicable
Expiry Date:	Not Applicable

Accessory Compatibility			
Wall Bracket	✓	Trolley Bracket	✓
Stand	✓	Sharpsafe Trolley**	✓
NPSD Tray	✗	Pedestal Trolley**	✓

Logistic Data	
Carton Quantity:	50
Carton Weight (when empty):	1.6kg
Carton Weight (when full):	23.5kg
Carton Dimensions (mm):	680 x 345 x 810 (HxWxD)
Carton Volume:	0.19m ³
Storage Temperature Conditions:	Ambient
Storage Lower Temp Limit:	None
Storage Upper Temp Limit:	None
Sterile/Non Sterile:	Non Sterile

Features and Benefits			
In-Mould Label	✓	Latex Free	✓
Temporary Closure	✓	Low Carbon Footprint Material	✓
Permanent Closure	✓	Leak Resistance	✓
Four Click Assembly	✓	Disconnection Features	✓
Fill Line	✓	Safety Shield	✓
Protected Access*	n/a	Semi Translucent Lid	✓
Single Use	✓	Carry Handle	✓
Different Lid Variations*	✓	Puncture Resistant	✓

Pallet Information	
Euro Pallet Carton Quantity:	9
Euro Pallet Dimensions (mm):	1200 x 800
Standard Pallet Carton Quantity:	12
Standard Pallet Dimensions (mm):	1200 x 1000

Manufacturing Information	
Made In:	Wales, UK
By:	Vernacare

*Weight without safety shield/Weight including safety shield
 **Different product codes may apply
 *** When used with a Safe Clip Trolley Bracket Adaptor