

RESCUE ESSENTIALS

# Compressed Gauze, 12ft, Z-Folded

Code: SR-GZ



- The unique folding and dispensing system makes the Compressed Z-Folded Gauze the ideal solution for fast and easy wound packing. The system allows this gauze to be easily used and cut with virtually no wasted gauze or time. In pre-hospital emergency care scenarios, rolled gauze can become damaged and lose its sterility.
- This compressed gauze maintains sterility in rugged, compact packaging that facilitates wound packing and wound dressing. In wound packing situations, gauze removal from the wound tract is made easier because this Compressed Gauze is a continuous-feed gauze. This compressed sterile fluff bandage is the perfect tool to use in conjunction with a compression bandage to deal with serious blood loss..

## Key Features:

- **12 feet, 6-ply cotton, 4.5" sterile** compressed gauze.
- **Vacuum sealed** for maximum portability.
- **Easy-open white tear notches** for quick and easy access for wound packing and wound dressing.
- **Unique z-fold dispensing system** resulting in controlled application with no wasted gauze.

