

## 4-WHEEL PATIENT FOLDING TRANSIT WHEELCHAIR

### FOR LEVEL GROUND

- ✓ MDD Compliant
- ✓ UKCA Marked & CE Marked



Primarily used to transport persons across level surfaces to a place of safety, our 4-Wheel Transit Chairs are similar to those used by hospitals and ambulance services and they can be used to carry patients either up or down the stairs. A minimum of two operators are required to carry the chair but additional may be necessary depending on the weights to be lifted. Because of the lifting and Manual Handling involved, it is recommended that full Risk Assessment and Personal Emergency Evacuation Plan (PEEP) should be carried out prior to purchase.

#### SPECIFICATION:

- Chair Weight 9.6kg
- Chair Length (open) 930mm
- Chair Length (folded) 830mm
- Chair Width 520mm
- Load capacity 120kg



#### AVAILABLE IN:



- ✓ Acts like a wheelchair for easy manoeuvrability when on level ground.
- ✓ The person feels secure once the safety belt is fastened.
- ✓ Flat seat, arms, and footrest help to relax the client.
- ✓ Helps to reduce the total weight when a person is in the chair.
- ✓ Handle and seat can be folded down to minimise the size.