

TRUECLOT®

Blood Simulant Concentrate

When diluted in water, TrueClot Blood Simulant Concentrate will create a highly realistic simulated blood for use in traumatic bleeding control, wound packing and advanced hemostatic dressings training. The simulated blood will form realistic simulated blood clots when used with TrueClot Simulated Hemostatic Gauze or TrueClot Clotting Solution.

Features

Small form factor, dilute with water when ready to use

Forms highly realistic, simulated blood clots when mixed with TrueClot Clotting Solution or Simulated Hemostatic Gauze

Washes off of skin with soap and water (for best washability use bar soap)

Washes out of clothing with household laundering

3 year shelf life

Available in 1 gallon (3.8 L) and 4 x 1 liter sizes



MADE IN THE USA

Dimensions

Gallon size: 6 x 9 x 1.5 in (15.2 x 22.3 x 3.8 cm) - 0.4 lbs (190 g)

4 x 1 liter size: 6 x 9 x 1.5 in (15.2 x 22.3 x 3.8 cm) - 0.5 lbs (225 g)

Item # 5115: 1 Gallon (3.8 L) packet - Item # 5117: 4 x 1 L packets

www.trueclot.com