

TRUECLOT®

CSAT Chest Seal Application Trainer

The TrueClot Chest Seal Application Trainer (CSAT) is both a table top and wearable scenario task trainer designed to teach chest seal application skills. Simulated wounds are used on the table top to demonstrate basic skills and worn on the vest to teach students to manage a patient while correctly placing a chest seal.

Features

The first task trainer that can be used as both a table top AND wearable prop

Simulated bubbling, breathing and chest seal venting/sealing

Wearable entry and exit wounds with active bleeding for scenario training

Durable simulated wounds and Velcro MOLLE vest

Apply and reuse real chest seals for realistic training

Easy to clean with water and air dry

Available in Dark and Light flesh tones



MADE IN THE USA

Dimensions

Kit in case: 11 x 13.5 x 5.2 in (28 x 34.3 x 13.2 cm) - 5.2 lbs (2.4 kg)

Wound surface: 6.25 x 8.75 in (15.9 x 22.2 cm) - Liquid capacity: 33.8 oz (1 liter)

Item # 5214

www.trueclot.com