

EXTREME SPORTS



RAPID  
APPLICATION

U.S. PATENTED  
TECHNOLOGY



POWERED BY





# DESIGNED TO WORK, **FAST.**

The RapidStop® Tourniquet is an award-winning, rapid application tourniquet that maximises the chances of survival. Developed at MIT and perfected by EMS, military and industry specialists. RapidStop® represents a haemostatic solution with military-grade ruggedness, all age suitability, superior single-handed application and a rapid application time.

Whether it's motorbike riding, mountain biking or rock climbing each of these activities occur in remote locations often inaccessible to emergency services. The unpredictable nature of the outdoors means disaster can strike anywhere at any time. Being prepared for the unexpected is crucial so it is important to have the right first aid kit for each situation.

RapidStop bleed control kits provide the necessary tools to save a life when disaster strikes. A rapid application tourniquet allows for quick response to major trauma injuries. With a trauma and haemorrhage control bandage as well as a RapidStop Tourniquet this kit provides the additional tools necessary to add to a standard workplace kit. Providing treatment options with efficiency and effectiveness in mind.



INTUITIVE

STRONG

LIFESAVING



## IF YOU NEED RapidStop®, YOU NEED IT FAST!

A severed brachial, radial or ulnar artery, presents a window of approximately 60 seconds before the decreasing blood pressure renders the injured unconscious. Shock can set in even sooner. Achieving full occlusion in the shortest time possible is critical. Faster application means faster stabilisation, increasing chances of survival.

## TO PLACE AN ORDER or MORE INFO:

[www.rapid-stop.com](http://www.rapid-stop.com)

[sales@aerohealthcare.com](mailto:sales@aerohealthcare.com)

AERO Healthcare | HQ

63 Seaton St, Armidale NSW 2350

1800 628 881



## INTUITIVE.

Fast and easy single-handed self-application are paramount to a tourniquet's effectiveness in any environment. RapidStop® stands out from the competition as one of the few tourniquets that exceed these life-saving requirements utilising key patented and proprietary technologies.

## STRONG.

Made in America and built tough to military standards. The patented ratchet lever technology is made from the same polycarbonate used in the F-22 jet fighter canopy. The tourniquet's band is made from high-strength nylon webbing rated to a two-ton breaking point. RapidStop® has constructed the tourniquet to endure extreme temperatures and resist any corrosion, ensuring it stands true in the toughest of conditions.

## LIFESAVING.

Easy as 1-2-3. Three quick and simple steps set RapidStop® Tourniquet apart from the competition. Clear instructions enable an untrained user to apply the tourniquet with speed and reliability. Integrated instructions on the RapidStop® Tourniquet are stitched on as weather-resistant waterproof vinyl tags. Pocket-sized instruction cards make learning and becoming familiar with the RapidStop® easy.